

Lifeline

MedReady® Medication Dispenser

Your patients rely on their medications to keep them healthy. But complicated medication schedules can lead to mistakes: missing doses, taking incorrect amounts, or taking pills at the wrong time. These mistakes could lead to unnecessary doctor or hospital visits.

We are pleased to offer a better way to help your patients manage their medications from home: **The MedReady medication dispenser.**

It's a simple, affordable and automated way to manage multiple medications with 24/7 monitoring.



Who would benefit from the MedReady medication dispenser?

MedReady is ideal for your patients who:

- Have recurring hospitalizations related to medication mistakes
- Have complicated medication regimens
- Take multiple medications
- Have a tendency to forget to take their medications

How does MedReady work?

1. **Fill:** Using the dosing templates (1 to 4 doses a day) add medications* from dose 1 then confirm when completed.
2. **Take medication:** The pill tray rotates to the scheduled dose and the audible alarm sounds. Open the medication compartment and retrieve the pills.
3. **Alerts:** Caregiver receives an SMS text or email¹ if medication has not been taken or when the dispenser needs troubleshooting.
4. **Lifeline Response Centre support:** Our specially trained Associates will assist the caregiver and can send the right help if needed.



* It is designed only for use with pills. Nitroglycerin, liquids, and powders should not be placed in the device.

Features and Benefits

Easy to use – Dispenses pre-loaded doses at the scheduled time². Simply press open the compartment to receive the medication.

Monitored by Lifeline's Response Centre – Added reassurance of knowing caregiver(s) are alerted when a dose is missed by SMS text or email. The Response Centre is also available to support the caregiver when needed.

Customizable – Can be programmed to dispense medication according to individual-specific prescriptions and holds up to 28 doses and accommodates up to 4 doses a day.

Highly reliable – Visual and audible alarms indicate when it is time to take the medication and allows for greater adherence.

Safe and secure – Tamper-proof construction, secure edges, and a safety lock help ensure that medications stay in the right hands.

Back up battery – Designed to distribute every dose on time, even in a power outage, thanks to our reliable backup battery³ protection.

For more information about MedReady or any of our Lifeline services, call:

1-800-LIFELINE
(1-800-543-3546)

MedReady® is a registered trademark of MedReady, Inc.

This service is only available at participating Lifeline programs.

1. The subscriber and/or caregiver require a smart phone in order to enroll in this service.

2. Certain limitations subject to third party cellular provider availability and coverage. Signal range may vary.

3. Back up power battery can last up to 72 hours.

www.lifeline.ca